

ANTHEM™

PLATINUM

Designed for vibrant lifestyles, includes 16 channels and 16 bands for optimal high-resolution sound imaging and frequency shaping — provides optimum flexibility and performance in a broad range of demanding listening environments that include high levels of background noise

WhistleFree Feedback Cancellation — the world's leading feedback cancellation system — enables open fit options on the widest range of hearing loss while providing intelligent artifact elimination

Premium **Environmental Noise Classification** optimizes 16-channel high resolution for maximum performance in noise as well as comfort control

Full **ClassiFi EA** includes maximum adjustability with 5 adaptation levels and precision with Speech, Speech & Noise, Machine Noise, Wind, Other Noise and Quiet environments

Industry's leading **directional system** features highest available directivity indexes, lowest operational noise floor and assured accuracy

Live Real Ear Measurement measures the hearing aid output in the patient's ear, ensuring an accurate fit to prescriptive targets

Live Speech Mapping verifies hearing aid's processing of speech, or any live acoustic input, in real time; 3D display engages patient and family in fitting process

Intuitive Features include **T¹**, Music Processing, Voice Indicators, Appointment Reminder and Automatic Telephone Solutions
Leisure Listening includes Music Processing with four genres and Television Processing

ANTHEM™

GOLD

Designed for active lifestyles, includes 12 channels and 12 bands for high-resolution sound imaging and frequency shaping — provides excellent performance in a broad range of demanding listening environments that include moderate levels of background noise

WhistleFree Feedback Cancellation — the world's leading feedback cancellation system — enables open fit options on a wide range of hearing loss while providing intelligent artifact elimination

Advanced Environmental Noise Classification utilizes 12-channel high resolution for successful performance in noise

ClassiFi EA includes adjustability with 4 adaptation levels and precision with Speech & Noise, Machine Noise, Wind, Other Noise and Quiet environments

Industry's leading **directional system** features highest available directivity indexes, lowest operational noise floor and assured accuracy

Live Real Ear Measurement measures the hearing aid output in the patient's ear, ensuring an accurate fit to prescriptive targets

Live Speech Mapping verifies hearing aid's processing of speech, or any live acoustic input, in real time

Intuitive Features include **T²** and Automatic Telephone Solutions
Leisure Listening includes Television Processing

ANTHEM: THE WIDEST RANGE OF ADVANCED HEARING SOLUTIONS AVAILABLE

ANTHEM™

SILVER

Designed for social lifestyles, includes 8 channels and 8 bands for precise sound imaging, providing good performance in environments with minimal to moderate levels of background noise

WhistleFree Feedback Cancellation — the world's leading feedback cancellation system — enables open fit options while providing intelligent artifact elimination

Select Environmental Noise Classification utilizes 8-channel resolution to improve performance in noise

ClassiFi EA is adjustable and includes Wind, Other Noise and Quiet environments
Industry's leading **directional system** features highest available directivity indexes

Live Real Ear Measurement measures the hearing aid output in the patient's ear

Live Speech Mapping verifies hearing aid's processing of speech, or any live acoustic input, in real time

Intuitive Features include **T¹** and Automatic Telephone Solutions
Leisure Listening includes Television Processing

ANTHEM™

BRONZE

Designed for calm lifestyles, includes 6 channels and bands for performance in environments with minimal to moderate levels of background noise

WhistleFree Feedback Cancellation — the world's leading feedback cancellation system
Standard Environmental Noise Classification utilizes 6-channel resolution for noise management

ClassiFi EA includes Wind, Other Noise and Quiet environments

Industry's leading **directional system** features highest available directivity indexes

Live Real Ear Measurement and **Live Speech Mapping**

Intuitive Features include **T¹**, Automatic Telephone Solutions and Television Processing